



COUCH CANE™



PRODUCT #2001



FITS MOST ANY COUCH, CHAIR, OR LIFT CHAIR




FEATURES



SIT AND STAND SAFELY




OPTIONAL TRAY
ACCESSORY ITEM #2055
sold separately



PROVIDES ASSISTANCE IN AND OUT OF COUCH OR CHAIR



SADDLE BAG POUCH



ADJUSTABLE IN HEIGHT AND DEPTH

SPECIFICATIONS

- Height Adjustment Handle: 34" to 40"
- Dimensions of Handle: 6" x 6" 1.25 diameter
- Base Length Adjustment: 20" to 36"
- Weight: 19.5 lbs
- Base Feet Measurements: 4" x 7"
- Weight Limit of Handle: 250 lbs
- Case of 3 Dimensions: 29" L x 19" W x 12" H; 60 lbs
- Package Dimensions: 29" L x 19" W x 4" H; 19 lbs