



# ASSIST-A-TRAY™



PRODUCT #2050



FITS MOST ANY COUCH, CHAIR, OR LIFT CHAIR




## FEATURES




SIT AND STAND SAFELY



INCLUDES CUP HOLDER AND UTENSIL COMPARTMENT



PROVIDES ASSISTANCE IN AND OUT OF COUCH OR CHAIR



MULTI-USE TRAY SWIVELS 360°



ADJUSTABLE IN HEIGHT AND DEPTH

## SPECIFICATIONS

- Height of Tray: 26" to 32"
- Height Adjustment Handle: 34" to 40"
- Dimensions of Handle: 6" x 6" 1.25 Diameter
- Base Length Adjustment: 20" to 36"
- Weight: 25 lbs
- Size of Tray: 16" x 18"
- Base Feet Measurements: 4" x 7"
- Weight Limit of Tray: 30 lbs
- Weight Limit of Handle: 250 lbs
- Case of 2 Dimensions: 29" L x 19" W x 8" H; 50 lbs
- Package Dimensions: 29" L x 19" W x 4" H; 25 lbs