



EZ STAND-N-GO™



PRODUCT #2200



FEATURES




SUPPORTS UP TO 300 LBS




FITS ON TOP OR UNDER CUSHION OF COUCH, CHAIR, OR RECLINER



UNIVERSAL - FITS MOST ANY COUCH, CHAIR, OR RECLINER



PROVIDES ASSISTANCE IN AND OUT OF COUCH OR CHAIR



ADJUSTABLE IN HEIGHT AND DEPTH



ERGONOMIC HANDLES FOR EXTRA COMFORT

SPECIFICATIONS

- **Chair Type:** Fully adjustable to fit most couches, chairs, and recliners
- **Handle Height From Floor:** Adjustable from 22" - 32"
- **Width Between Handles:** Adjustable from 19" - 26"
- **Base Dimensions:** Adjustable Depth 23" - 26" & Adjustable Width 19" - 26"
- **Tube Diameter:** 1"
- **Product Dimensions:** 23" - 26" x 19" - 26" x 22" - 32" 8.5 lbs
- **Package Dimensions:** 21" x 12.5" x 5.5" 8.5 lbs
- **Case of 5 Dimensions:** 5.5" x 21.25" x 12.5" 42.5 lbs