



OMNI TRAY™



PRODUCT #2600



FITS MOST ANY COUCH, CHAIR, OR LIFT CHAIR



ACCESSORIES (SOLD SEPARATELY)



UTENSIL TRAY
PRODUCT # 2602



CUP HOLDER
PRODUCT # 2601



ORGANIZER
PRODUCT # 2603

FEATURES




LARGE, OVERSIZED TRAY



TRAY PIVOTS OUT OF THE WAY WHEN NOT IN USE




SIT AND STAND SAFELY



PROVIDES ASSISTANCE IN AND OUT OF COUCH OR CHAIR



MULTI-USE TRAY SWIVELS 360°



ADJUSTABLE IN HEIGHT AND DEPTH

SPECIFICATIONS

- Height of Tray: 24" – 32"
- Height of Handle Adjustment: 28.5" x 36.25"
- Dimensions of Handle: 4.5" x 5"
- Base Length Adjustment: 25.5" - 35.5"
- Weight: 20 lbs
- Size of Tray: 20" x 15"
- Base Feet Measurements: 8" x 5"
- Weight Limit of Tray: 30 lbs
- Weight Limit of Handle: 250 lbs
- Case of 2 Dimensions: 20.5" x 21" x 7" – 40 lbs
- Package Dimensions: 20.5" x 21" x 3.5" – 20 lbs