



EZ STAND-N-GO HD™



PRODUCT #9220




FEATURES




SUPPORTS UP TO 400 LBS



FITS ON TOP OR UNDER CUSHION OF COUCH OR CHAIR



PROVIDES ASSISTANCE IN AND OUT OF COUCH OR CHAIR



ADJUSTABLE IN HEIGHT AND DEPTH



ERGONOMIC HANDLES FOR EXTRA COMFORT

SPECIFICATIONS

- **Chair Type:** Fully adjustable to fit most couches, chairs, and recliners
- **Handle Height From Floor:** Adjustable from 22" - 32"
- **Width Between Handles:** Adjustable from 19" - 26"
- **Base Dimensions:** Adjustable Depth 23" - 26" & Adjustable Width 19" - 26"
- **Tube Diameter:** 1"
- **Product Dimensions:** 23" - 26" x 19" - 26" x 22" - 32"; 8.5 lbs
- **Package Dimensions:** 21" x 12.5" x 5.5" 9.5 lbs
- **Case of 5 Dimensions:** 25.5" x 21.25" x 12.5"; 48.5 lbs