



BEDRAIL ADVANTAGE TRAVELER™



PRODUCT #5000

FEATURES



INCLUDES ORGANIZER POUCH




EXCELLENT FOR BARIATRIC USE
SUPPORTS UP TO 400 LBS



FOLDS FOR TRAVEL
OR STORAGE



PROVIDES ASSISTANCE
IN AND OUT OF BED



ONLY USE WITH 10" TO 16" MATTRESSES; NOT
INTENDED FOR USE ON ADJUSTABLE BEDS
OR BEDS WITH EXTRA SOFT MATTRESSES



SECURES TO BED FRAME WITH
INCLUDED SAFETY STRAP

SPECIFICATIONS

- Rail Height: 22"
- Rail Width: 17"
- Base Dimensions: 21" x 17"
- Weight of Product: 8 lbs
- Weight Capacity: 400 lbs
- Package Dimensions: 23" x 18" x 2"
- Case of 4 Dimensions: 23" x 18" x 8"; 34 lbs