

BEDSIDE ECONORAIL™







SNAP TOGETHER ASSEMBLY




WEIGHS LESS THAN 2 LBS, SUPPORTS UP TO 300 LBS



INCLUDES ORGANIZER POUCH



PROVIDES ASSISTANCE IN AND OUT OF BED



ONLY USE WITH 8" TO 16" MATTRESSES; NOT INTENDED FOR USE ON ADJUSTABLE BEDS OR BEDS WITH EXTRA SOFT MATTRESSES



SECURES TO BED FRAME WIT INCLUDED SAFETY STRAP

- Handle Height: 21.5"
- Handle Width: 16" at base to 12" at top
- Base Dimensions: 20" Deep × 16" Wide
- Tube Diameter: I"
- Weight of Product: 1.5 lbs
- Made of: Anodized aluminum
- Package Dimensions: 17.5" L × 12" W × 1.625" H
- Case of 5 Dimensions: 17.5" L × 12" W × 5" H; 6 lbs

PRODUCT #5100