



STABLE RAIL™



PRODUCT #5800



FEATURES



ADJUSTABLE LEGS
PROVIDE STABILITY



SUPPORTS UP TO 300 LBS




INCLUDES ORGANIZER POUCH



ACCESSORIZE THE STABLE RAIL
WITH THE STANDER ARM #5810



PROVIDES ASSISTANCE
IN AND OUT OF BED



HEIGHT ADJUSTS TO FIT 8"-18" MATTRESSES;
NOT INTENDED FOR USE ON ADJUSTABLE BEDS
OR BEDS WITH EXTRA SOFT MATTRESSES



SECURES TO BED FRAME WITH
INCLUDED SAFETY STRAP

SPECIFICATIONS

- Height Adjustment From Floor to Top of Handle: 21"- 44"
- Rail Width: 19"
- Height Range From the Floor to Horizontal Support Structure Under Mattress: 13.5"- 21.5"
- Weight of Product: 10.5 lbs
- Package Dimensions:
28" L x 20.25" W x 2" H
- Case of 4 Dimensions:
28" L x 20.25" W x 5.75" H; 32 lbs