



# 30" SAFETY BED RAIL™



PRODUCT #8050



## FEATURES



SAFETY BED RAIL PIVOTS OUT OF THE WAY



SUPPORTS UP TO 300 LBS



OPTIONAL PADDED POUCH  
AVAILABLE FOR  
PURCHASE #8053



PROVIDES ASSISTANCE  
IN AND OUT OF BED



ONLY USE WITH 10" TO 19" MATTRESSES; NOT  
INTENDED FOR USE ON ADJUSTABLE BEDS OR  
BEDS WITH EXTRA SOFT MATTRESSES



SECURES TO BED FRAME WITH  
INCLUDED SAFETY STRAP

## SPECIFICATIONS

- Width of Rail: 30"
- Height of Rail: 23" high from the base of mattress
- Gap Between Rails: 3"
- Weight of Product: 12 lbs
- Package Dimensions: 32" L x 21.5" W x 1.25" H 12 lbs
- Case Dimensions: 32" L x 21.5" W x 5" H; 50 lbs