



PRIME SAFETY BED RAIL™



PRODUCT #8940



FEATURES



INCLUDES ORGANIZER POUCH



HOOKS TO HANG ON FOOT OR SIDE OF BED



SUPPORTS UP TO 400 LBS.



ATTACH PRODUCT #8905 TO EXTEND LENGTH OF RAIL



PROVIDES ASSISTANCE IN AND OUT OF BED



ONLY USE WITH 6" TO 15" MATTRESSES



SECURES TO BED FRAME WITH INCLUDED SAFETY STRAP

SPECIFICATIONS

- Rail Height: 18.5"
- Rail Width: 27"
- Base Dimensions: 19.5" x 28.5"
- Weight of Product: 15.7 lbs
- Package Dimensions: 21.25" x 32" x 4.25"
- Case of 2 Dimensions: 22" x 32.5" x 9"; 35.10 lbs