



NEW

EZ CLICK LTC BED HANDLE™

HEAVY DUTY
UP TO 500 LBS.



PRODUCT #8950



OVERCOMES
ALL 7 FDA
ENTRAPMENT ZONES

FEATURES



ATTACHES TO MOST
INSTITUTIONAL BEDS



SUPPORTS UP TO 500 LBS



2-SIDED BED HANDLE



DESIGNED TO PREVENT ENTRAPMENT



EZ CLICK BUTTON RELEASE



PROVIDES ASSISTANCE
IN AND OUT OF BED



ONLY USE WITH 6" TO 9" MATTRESSES;
NOT INTENDED FOR BEDS WITH EXTRA SOFT MATTRESSES



ANTI-MICROBIAL

SPECIFICATIONS

- *Width of Rail: 20.5" after installed on bed*
- *Height of Rail: 13" high from the base of mattress*
- *Weight of Product: 18.39 lbs*
- *Package Dimensions: 27" x 13" x 4.72"*
- *Case of 2 Dimensions: 27.5 L" X 13.75" W X 10.5" H; 46.75 lbs*