



BATHTUB SECURITY POLE & CURVE GRAB BAR™

NEW



FEATURES



CLAMPS TO 3" – 7" EDGE OF TUB




SUPPORTS UP TO 300 LBS



TENSION MOUNTED – NO SCREWS OR DRILLING



ADJUSTS TO FIT 7' – 9' CEILINGS



MAKES SITTING AND STANDING EASY



LADDER GRIP ROTATES & LOCKS EVERY 45°

SPECIFICATIONS

- **Materials:** Rust resistant zinc plated steel with durable powder coat finish
- **Fits Ceiling Heights:** Range of 7'-9' ft.
- **Weight of Product:** 22 lbs
- **Package Dimensions:** 55.9" L x 13.5" W x 2.95" H; 24 lbs
- **Case of 2 Dimensions:** 56.69" L x 14.17" W x 7.28" H; 52 lbs

COLOR

WHITE #1120-W

PRODUCT #1120

EZ SWIVEL OVERBED TABLE™

NEW



FEATURES



TRAY ROTATES 360° WHEN NOT IN USE



INCLUDES ORGANIZER POUCH



LARGE, OVERSIZED TRAY



INCLUDES CUP HOLDER



ADJUSTS IN HEIGHT TO FIT MOST ANY BED

SPECIFICATIONS

- **Height of Tray:** 24" – 32"
- **Height of Handle Adjustment:** 28.5" x 36.25"
- **Base Length Adjustment:** 25.5" - 35.5"
- **Weight:** 20 lbs
- **Size of Tray:** 20" x 15"
- **Base Measurements:** 6"
- **Weight Limit of Tray:** 30 lbs
- **Case of 2 Dimensions:** 20.5" x 21" x 7" – 40 lbs
- **Package Dimensions:** 20.5" x 21" x 3.5" – 20 lbs

PRODUCT #6000