



EZ FOLD-N-GO WALKER - LITE™

EZ FOLD-N-GO WALKER - LITE+™

COLOR



Black Walnut
#4370-BW

FEATURES



INCLUDES ORGANIZER POUCH



RUBBER FEET ON LEGS



WEIGHS ONLY 6 LBS
SUPPORTS UP TO 300 LBS



HANDLE HEIGHT ADJUSTMENT:
32.5" - 39"



PRESS TO LOCK
LIFT TO UNLOCK

PRODUCT #4370

SPECIFICATIONS

- Weight: 6 lbs
- Collapsed Dimension: 36" H x 5.25" W x 8" D
- Open Depth: 21.5"
- Open Width: 25.5"
- Open Handle Width: 18"
- Package Dimensions: 32.6"L x 6.69"W x 4.72"H, 8.2 lbs
- Case Dimensions: 33.5"L x 7.28"W x 15"H, 24.6 lbs
- Recommended body height range: 5'6" - 6'3"

COLOR



Black Walnut
#4380-BW

FEATURES



INCLUDES ORGANIZER POUCH



5" WHEELS AND SKI GLIDES



WEIGHS ONLY 7 LBS
SUPPORTS UP TO 300 LBS



HANDLE HEIGHT ADJUSTMENT:
32.5" - 39"



PRESS TO LOCK
LIFT TO UNLOCK

PRODUCT #4380

SPECIFICATIONS

- Weight: 7 lbs
- Collapsed Dimension: 36" H x 5.25" W x 8" D
- Wheels: 5" diameter wheels
- Case Dimensions: 37.25" L x 9.5" W x 19.25" H, 28 lbs
- Open Depth: 21.5"
- Open Width: 25.5"
- Open Handle Width: 18"
- Package Dimensions: 36.25"L x 8.5"W x 6"H, 8 lbs
- Case Dimensions: 37.25"L x 9.5"W x 19.25"H, 28 lbs
- Recommended body height range: 5'6" - 6'3"